




## Introducing the V.A.C. VERAFLU CLEANSE CHOICE™ Dressing for V.A.C. VERAFLU™ Therapy

A WOUND CLEANSING OPTION FOR CLINICIANS WHEN SURGICAL DEBRIDEMENT MUST BE DELAYED OR IS NOT POSSIBLE OR APPROPRIATE




	Innovative Features	Benefits
	<p><i>Unique dressing:</i></p> <ul style="list-style-type: none"> <li>• Provides "mechanical" movement at the wound surface in combination with cyclic delivery and dwell of topical solutions</li> <li>• Facilitates removal of thick wound exudate, such as fibrin, slough and other infectious materials</li> </ul>	<ul style="list-style-type: none"> <li>• Can be used early to initiate wound cleansing therapy</li> <li>• Ideal for dirty wounds needing immediate wound cleansing</li> </ul>
	<p><i>3-layer foam design includes:</i></p> <ul style="list-style-type: none"> <li>• A contact layer with a pattern of 1 cm through holes</li> <li>• Two cover layers (without holes) to provide application options for wounds of varying depths</li> </ul>	<ul style="list-style-type: none"> <li>• Multiple layers provide improved conformability</li> <li>• Cover layer is designed to capture large wound exudate particles before they enter into the tubing</li> <li>• Single pad application</li> </ul>

## V.A.C. VERAFLU CLEANSE CHOICE™ DRESSING: A CONVENIENT PATH TO WOUND CLOSURE

### PRELIMINARY EVIDENCE SUGGESTS:<sup>1</sup>

The adjunctive use of V.A.C. VERAFLU CLEANSE CHOICE™ Dressing for V.A.C. VERAFLU™ Therapy may be suitable for wound cleansing in wounds when complete surgical debridement is not possible or appropriate and areas of non-viable tissue are still present on the wound surface.

- Assisted in loosening, solubilizing, detaching and removing viscous exudate, dry fibrin, wet slough and other infectious materials in a retrospective review of the outcomes of 21 patients with large complex chronic wounds
- Percent surface area of black non-viable tissue and yellow fibrinous slough was reduced to ≤10% in 18/21 (85.7%) and 12/21 (57.1%) of wounds after an average of 1-3 applications (3-9 days) using saline solution



1. Photos courtesy of Teot L., Boissiere F., Fluieraru, S. Novel foam dressing using negative pressure wound therapy with instillation to remove thick exudate. *International Wound Journal* 2017; 10:1111/iwj.12719.



**New V.A.C. VERAFLU CLEANSE CHOICE™ Dressing MOA Video** highlights the V.A.C. VERAFLU™ Therapy difference. See it today!

#### Ordering Information

SKU	Description	Size	Quantity
ULTVCC05MD	V.A.C. VERAFLU CLEANSE CHOICE™ Dressing	Medium	5-Pack

For more information, contact your local Acelity™ Representative.

**NOTE: Specific indications, contraindications, warnings, precautions, and safety information exist for KCI products and therapies. Please consult a physician and product instructions for use prior to application. This material is intended for healthcare professionals only.**