



A SIMPLE SOLUTION  
FOR COMPLEX WOUNDS

INTRODUCING

# V.A.C. VERAFLU CLEANSE CHOICE™ Dressing for V.A.C. VERAFLU™ Therapy




- The V.A.C. VERAFLU CLEANSE CHOICE™ Dressing, in conjunction with V.A.C. VERAFLU™ Therapy, can be used to initiate therapy and facilitate removal of infectious materials.
- V.A.C. VERAFLU CLEANSE CHOICE™ provides a wound cleansing option for clinicians when surgical debridement must be delayed or is not possible or appropriate.<sup>1</sup>

# How does the V.A.C. VERAFLU CLEANSE CHOICE™ Dressing work with V.A.C. VERAFLU™ Therapy?

- Dressing includes 3 different layers: a contact layer with a pattern of larger through holes and 2 cover layers to provide application options for wounds with varying depths
- Provides “mechanical” movement at the wound surface in combination with cyclic delivery and dwell of topical solutions<sup>2</sup>
- The contact layer in combination with a cover layer captures larger particles to help prevent blockage alarms
- The contact layer and the two cover layers provide enhanced conformability

The V.A.C. VERAFLU CLEANSE CHOICE™ Dressing has been designed to help facilitate the removal of infectious material such as thick fibrinous exudate, slough and other wound bioburden up to 30 centipoise(cP).<sup>3</sup>



These are general averages for reference only.



# When to use V.A.C. VERAFL0™ Therapy



Physician panel recommends that when used in conjunction with good clinical practice, V.A.C. VERAFL0™ Therapy may be used in complex wounds such as:

- Contaminated wounds
- Diabetic wounds
- Acute and chronically infected wounds
- Traumatic wounds
- Pressure ulcers
- As adjunctive therapy in infected wounds in the presence of orthopaedic hardware or joint implants

## CHOOSING V.A.C. VERAFL0™ THERAPY IS AS EASY AS:

- 1** Choose your V.A.C. VERAFL0™ Dressing (see back page)
- 2** Choose your settings (dwell time and cycle frequency)
- 3** Choose your solution\*

Physician panel recommendations based on published evidence and their experience:<sup>4</sup>



**DWELL TIME**  
10-20 Minutes



**CYCLE FREQUENCY**  
2-4 Hours at -125mmHg  
(up to 6 hours for larger wounds)



**SOLUTION**  
Saline  
(in most instances)

NPWTi-d is an adjunct to other standard principles of appropriate wound assessment and treatment.

## ORDERING INFORMATION:

### Description:

V.A.C. VERAFL0  
CLEANSE CHOICE™ Dressing

### SKU:

ULTVCC05MD  
ULTVCC05LG

### Quantity:

5-Pack

\*Please contact your KCI representative for a list of topical wound solutions compatible with KCI dressings and components.





	V.A.C. VERAFLOR™ Dressings	V.A.C. VERAFLOR™ Large Dressing	V.A.C. VERAFLOR CLEANSE™ Dressing	V.A.C. VERAFLOR CLEANSE CHOICE™ Dressing
--	----------------------------	---------------------------------	-----------------------------------	--

	V.A.C. VERAFLOR™ Dressings	V.A.C. VERAFLOR™ Large Dressing	V.A.C. VERAFLOR CLEANSE™ Dressing	V.A.C. VERAFLOR CLEANSE CHOICE™ Dressing
<b>Wound Characteristics</b>	Open wounds, including wounds with shallow undermining or tunnel areas where the distal aspect is visible	Large open wounds, including wounds with shallow undermining or tunnel areas where the distal aspect is visible	Cavity wounds or wounds with complex geometries	Wounds with thick fibrinous exudate, slough, infectious material and other wound bioburden
<b>Key Goal(s) of Therapy</b>	<ul style="list-style-type: none"> <li>When used in conjunction with V.A.C. VERAFLOR™ Therapy, to help facilitate the removal of infectious material and other wound bioburden.</li> <li>Generation of robust granulation tissue</li> </ul>	<ul style="list-style-type: none"> <li>When used in conjunction with V.A.C. VERAFLOR™ Therapy, to help facilitate the removal of infectious material and other wound bioburden.</li> <li>Generation of robust granulation tissue in large wounds</li> </ul>	<ul style="list-style-type: none"> <li>When used in conjunction with V.A.C. VERAFLOR™ Therapy, to initiate therapy and to help facilitate the removal of infectious material and other wound bioburden.</li> <li>Easy application into explored tunnels and undermining</li> </ul>	<ul style="list-style-type: none"> <li>When used in conjunction with V.A.C. VERAFLOR™ Therapy, to initiate therapy and to help facilitate the removal of infectious material such as thick fibrinous exudate, slough, and other wound bioburden.</li> <li>Provides a wound cleansing option for clinicians when surgical debridement must be delayed or is not possible or appropriate.<sup>1</sup></li> </ul>
<b>Shape</b>	Spiral-cut foam	Block foam pre-slit into two layers	Tubular shape	Block foam pre-slit into three layers
<b>Application Characteristics</b>	Easy application: <ul style="list-style-type: none"> <li>Size without scissors</li> <li>Precut area for pad application when used for bridging</li> <li>Single pad application</li> </ul>	Easy application: <ul style="list-style-type: none"> <li>Ideal for large surface areas with shallow depths</li> <li>Provided with V.A.C. VERAT.R.A.C. DUO™ Tube set for extended surface area coverage</li> </ul>	Application flexibility: <ul style="list-style-type: none"> <li>Ideal for addressing wounds with complex geometries (e.g. tunnels, undermining and circumferential applications)</li> <li>Single pad application</li> </ul>	Easy application with flexibility: <ul style="list-style-type: none"> <li>Thin layers for improved conformability</li> <li>Multiple layers provide application options for wounds with varying depths</li> <li>Ideal for dirty wounds needing active therapy</li> <li>Single or duo pad application</li> </ul>
<b>Dressing Characteristics</b>	<p><b>Small:</b> 7.7 x 11.2 x 1.75cm per piece (2 pieces spiral)</p> <ul style="list-style-type: none"> <li>Each piece 24cm uncoiled</li> <li>Each piece 18cm long bridge</li> </ul> <p>- 2x V.A.C.® Advanced Drape - 1x V.A.C. VERAT.R.A.C.™ Pad - 2x 3M™ Cavilon™ No Sting Barrier Film</p> <p><b>Medium:</b> 14.7 x 17.4 x 1.75cm per piece (2 pieces spiral)</p> <ul style="list-style-type: none"> <li>Each piece 80cm uncoiled</li> <li>Each piece 74cm long bridge</li> </ul> <p>- 2x V.A.C.® Advanced Drape - 1x V.A.C. VERAT.R.A.C.™ Pad - 2x 3M™ Cavilon™ No Sting Barrier Film</p> <p><b>FILL ASSIST SET POINT RANGE</b></p> <p><b>Small:</b> 12-80mL (1 piece) 26-160mL (2 pieces)</p> <p><b>Medium:</b> 38-240mL (1 piece); 80-450mL (2 pieces)</p>	<p>25.0 x 15.0 x 1.6 cm per piece (oval - 2 pieces)</p> <p>- 5x V.A.C.® Advanced Drape - 1x V.A.C. V.A.C. VERAT.R.A.C. Duo™ Pad - 4x 3M™ Cavilon™ No Sting Barrier Film</p> <p><b>FILL ASSIST SET POINT RANGE</b></p> <p>55-350mL (1 piece); 110-500mL (2 pieces)</p>	<p>1 each coiled rod (24" rod; 1.25" diameter)</p> <p>- 2x V.A.C.® Advanced Drape - 1x V.A.C. VERAT.R.A.C.™ Pad - 2x 3M™ Cavilon™ No Sting Barrier Film</p> <p><b>FILL ASSIST SET POINT RANGE</b></p> <p>65-130mL (1/2 rod); 130-250mL (full rod)</p>	<p><b>Medium:</b> 18.0cm x 12.5cm x Varying Thicknesses (0.8cm – 3.2cm)</p> <p>- 4x V.A.C.® Advanced Drape - 1x V.A.C. VERAT.R.A.C.™ Pad - 4x 3M™ Cavilon™ No Sting Barrier Film</p> <p><b>Large:</b> 25.6cm x 15.0cm x Varying Thicknesses (0.8cm - 3.2cm)</p> <p>- 6x V.A.C.® Advanced Drape - 1x V.A.C. VERAT.R.A.C.™ Duo Pad - 5x 3M™ Cavilon™ No Sting Barrier Film</p> <p><b>FILL VOLUME START POINTS</b></p> <p><b>Medium:</b> 85mL (1.6cm Cover Layer); 42mL (0.8cm Cover Layer); 24mL (0.8cm Wound Contact Layer)</p> <p><b>Large:</b> 150mL (1.6cm Cover Layer); 75mL (0.8cm Cover Layer); 42mL (0.8cm Wound Contact Layer)</p>

**For more information, contact your local KCI Representative**

1. Teot L, Boissiere F, Fluieraru S. Novel foam dressing using negative pressure wound therapy with instillation to remove thick exudate. Int Wound J. 2017; doi: 10.1111/iwj.1217.  
 2. KCI. V.A.C. VERAFLOR™ Cleanse Choice™ Dressing Disruption of Simulated Tissue Test. DOC-0000039656. June 2015  
 3. KCI. V.A.C. VERAFLOR CLEANSE CHOICE™ Thick Exudate Fluid Removal & Pressure Maintenance Test Report. DOC-0000039947. June 2015.  
 4. Kim PJ, et al. Negative Pressure Wound Therapy With Instillation: Review of Evidence and Recommendations. Wounds. 2015 Dec;27(12):S2-S19.

**NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for KCI products and therapies. Please consult a physician and product instructions for use prior to application. This material is intended for healthcare professionals.**

©2016, 2018 KCI Licensing, Inc. 3M and Cavilon No Sting Barrier Film are trademarks of 3M Company. Unless otherwise designated, all trademarks are proprietary to KCI Licensing, Inc., its affiliates and/or licensors. All rights reserved. PRA-PM-AT-00023 (11/18)

