

PROMOGRAN™

MATRIX WOUND DRESSING



What is it?

PROMOGRAN™ Matrix is a sterile primary dressing made of freeze-dried composite of 45% oxidized regenerated cellulose (ORC) and 55% collagen. In the presence of exudate, the PROMOGRAN™ Matrix transforms into a soft, conformable, biodegradable gel, and thus allows contact with all areas of the wound.

What is it for?


PROMOGRAN™ Matrix, when covered with a semi-occlusive dressing, like a foam, maintains a moist wound environment. This environment is conducive to granulation tissue formation, epithelialization and optimal wound healing.

How to apply

Preparing the wound


Prepare the wound according to appropriate wound management protocol as instructed by your clinician. Wash your hands before and after touching the wound. Wear gloves if you have them.

Dressing preparation




1. Peel back the tray cover and leave dressing in tray.
2. If necessary, PROMOGRAN™ Matrix can be cut using clean scissors, ripped or folded to fit the wound bed.
3. For a small or deep exuding wound, cut the dressing prior to application or pack dressing into the wound.

Dressing application



1. If the wound bed is wet apply the dressing directly.
2. If the wound is dry moisten the dressing with saline (use the dressing tray as a dish to pre-wet)

Secondary dressing application



Cover with an appropriate secondary dressing as directed by your clinician.

How often to apply

It is not necessary to remove any remaining PROMOGRAN™ Matrix dressing during secondary dressing changes. Apply a new PROMOGRAN™ Matrix as directed by your clinician.

Frequently Asked Questions

Do I have to remove PROMOGRAN™ Matrix?

No, it is not necessary to remove PROMOGRAN™ Matrix. You do however need to remove the secondary dressing as directed by your clinician.

NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. The instructions for use are intended for healthcare professionals.

NOTE: This document is intended for patient self-care as guided by their HCP and is not intended to replace the Instructions for Use.

©Copyright 2020 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited. PRA-PM-ALL-00363 (04/20)

