

ALL YOU COULD WISH FOR IN A SILICONE CONTACT LAYER

WHAT IS IT?

ADAPTIC TOUCH® is a flexible, open-mesh silicone wound contact layer designed to optimise fluid management and minimise pain at dressing change.

COMFORT

ADAPTIC TOUCH® helps protect the wound while minimising patient pain and trauma at dressing change¹, raising both patients and clinicians' satisfaction².



PERFORMANCE

ADAPTIC TOUCH® open mesh structure allows free passage of exudate into the absorbent secondary dressing, reducing the risk of exudate pooling and maceration²⁻⁴.




The small size of the pores minimises the risk of granulation tissue in-growth⁴.

AFFORDABILITY

Treatment with ADAPTIC TOUCH® has been shown to be cost effective in clinical practice⁵.



ADAPTIC TOUCH® is up to 29% cheaper than Mepitel*.



ADAPTIC TOUCH®

NON-ADHERING SILICONE DRESSING




USER GUIDE

BEFORE APPLICATION


Refer to the instructions for use. Select a size of ADAPTIC TOUCH® that is larger than the wound. Prepare the wound according to wound management protocol. Ensure skin surrounding the wound is dry. Remove ADAPTIC TOUCH® from the pouch.

Dressing preparation



1. If needed **ADAPTIC TOUCH®** may be cut to size with sterile scissors.
2. It is best to cut **ADAPTIC TOUCH®** whilst one or two backing papers are in place.

Dressing application



1. Remove one backing paper from **ADAPTIC TOUCH®**. Place **ADAPTIC TOUCH®** dressing directly over the wound and smooth in place around the wound.
2. Remove the second backing paper using surgical tongs or forceps. If more than one piece of **ADAPTIC TOUCH®** is required, ensure dressings overlap to avoid secondary dressing adherence to the wound. Overlap should be minimised to prevent occlusion of holes.

Secondary dressing application



Cover **ADAPTIC TOUCH®** with a secondary dressing such as **TIELLE®**.

Choose the secondary dressing based on the level of exudate. Apply any secondary cover or retention/compression bandage.

Dressing change and removal



Dressing change frequency is dictated by good wound care practice and will depend on the condition of the wound. **ADAPTIC TOUCH®** may be left in place for several days.

1. Gently remove the secondary dressing.
2. Gently remove the dressing from the wound bed and discard.

WHEN TO USE

Management of dry to heavily exuding, partial and full-thickness chronic wounds, including:

- Leg ulcers (e.g. venous, arterial, mixed)
- Pressure ulcers
- Diabetic ulcers
- Donor sites
- 1st and 2nd degree burns
- Traumatic (e. g. skin tears)
- Surgical wounds



Suitable for use, under medical supervision, with negative pressure wound therapy (NPWT)



ADAPTIC TOUCH® can be used in patients with fragile or friable skin²

In vivo study comparing ADAPTIC TOUCH® performance to Mepitel showed that...

The smaller pore size of ADAPTIC TOUCH® may mean reduced risk of the secondary dressing adhering to the wound bed through the primary layer, which may lead to tissue damage and pain during dressing change and removal⁴.

ADAPTIC TOUCH®

Size	Dressings per box	Systagenix code
5cm x 7.6cm	10	TCH501
7.6cm x 11cm	10	TCH502
12.7cm x 15cm	10	TCH503
20cm x 32cm	5	TCH504

To learn more about the benefits of ADAPTIC TOUCH® contact your Systagenix representative at +44 (0) 20 302 787 16 or visit www.systagenix.com

References

* Based on UK drug tariff February 2014.

1. Ivin, N., and Harding, K. G. A series of case studies using a non-adhering silicone dressing. Poster, EWMA 2011. 2. International case series: Using ADAPTIC TOUCH® Non-Adhering Silicone Dressing: Case Studies. London: Wounds International, 2013. 3. Stephens S. et al. Evaluation of a non-adhering silicone wound contact dressing with optimised design for the management of dry to heavily exuding wounds. Poster, WIC 2011. 4. Brosnan, P. et al. The design principles of non-adherent materials. Poster, EWMA 2012. 5. Bianchi, J. Consensus guidance for the use of ADAPTIC TOUCH® non-adherent dressing. Wounds UK 2011.

© Systagenix wound management 2015. Brands marked with ® or ™ are trademarks of Systagenix. All other products referenced herein are acknowledged to be trademarks of their respective owners.